

## ● SUCTION TABLE TIMOBIL

The mobile suction station a big plus in equipment



### Technical information

Variable, well thought-out and especially strong in suction: With the wear-and-tear- and maintenance-free Timobil you have at your disposal a mobile suction unit which leaves absolutely nothing to be desired, whether in the surgery or in the laboratory. Naturally the Timobil – like almost all Wassermann suction units – works according to the low pressure principle. That means: The suction takes place at a very low air speed. The advantage for you: it is relaxed to work with, because of its virtual noiselessness. A further plus point is the high suction volume whereby, even in continuous use, the hands will not become cold. And by virtue of the automatic on switch the Timobil can be especially practically started via the ergonomic and vibration-free motor handpiece.

- mobile, maintenance- and wear-and-tear-free suction unit for the surgery and the laboratory
- with robust and ergonomic equipment engine
- strong torque, high rotation accuracy, brushless motor
- runs quietly, suitable for continuous operation
- no heat eduction, no cooling down of the hands
- perspex protection disc in front of the work area
- high protection via time-delay automatics of the suction unit
- high quality, quickly changed filterbag

## ● Technical data and accessories

Timobil	Item no.: 149999
Voltage*	220–240 V / 50/60 Hz
Power consumption	1.5 A
Female connector power consumption max	8 A
Output	350 W
W x H x D	530 x 900 x 480 mm
Weight	28.1 kg
Suction capacity	500 m <sup>3</sup> / h (free-blowing)
Control unit capacity	160 W
Handpiece motor capacity	95 W
Speed handpiece motor	1000 – 40 000 min <sup>-1</sup>
Duty	S6 (4 min. / 10 min.)

\*Other voltages on request

Incl. in delivery	Item no.:
Motor handpiece	111121
Rest for handpiece	111103
Protection sheet	149035
Micro fine dust filterbag	---
Accessories	Item no.:
Extractor hood	149130
Extractor hood with illumination	149140
Spare filterbag for micro fine dust (dust class M), 1 set = 4 pieces	611010
Spare filterbag for plaster (class 1 / class 2), 1 set = 4 pieces	611018
Storage tray	149010

Technical changes reserved



Timobil with special accessory extractor hood with illumination



Timobil with special accessory extractor hood



Motor handpiece

### Major features at a glance

- mobile suction unit for the surgery and the laboratory
- ergonomic and high-performance motor handpiece
- virtually noiseless operation, no cooling down of the hands
- suction unit can be started via motor handpiece
- completely free of maintenance and wear-and-tear

[www.wassermann-dental.com](http://www.wassermann-dental.com)