

## ● STEAM-JET WASI-STEAM 2

Absolutly clean



Wasi-Steam 2 with accessory splash protection

### Technical information

Strong, stronger, Wasi-Steam 2: This sophisticated steam-jet is the most powerful unit on the market. Thorough and efficient, it is suitable for all cleaning in the lab and in the field. Its high cleaning pressure of up to 6 bar also guarantees a high cleaning temperature. The steam quantity and wet steam saturation are easily regulated using the swivel head. The permanent water connection makes manual filling unnecessary. And it operates without chemical additives and removes limescale electronically, so it protects the environment.

- powerful steam-jet for all cleaning applications
- user-friendly operation using swivel head to regulate wet steam saturation and steam quantity
- steam pressure 5 bar at rest / 6 bar in operation
- protects seals by reducing pressure to 5 bar after 10-minute break
- selector for wet or dry steam
- water connection, electronic delimer, service interval display and huge stainless steel boiler

## ● Technical data and accessories

Wasi-Steam 2	Item no.: 174992	Incl. in delivery	Item no.:
Voltage*	220-240 V / 50/60 Hz	Discharge hose	320031
Power consumption	8.3 A (Si 10 A)	Supply hose	320060
Output	1880 W	Square socket spanner	174133
W x H x D	450 x 490 x 265 mm	<b>Accessories</b>	
Weight	30 kg	Splash protection (W max. 555 x H max. 250 x D max. 160 (left) / 210 (right) mm)	174230
Sound level	≤ 70 dB (A)	Fastening set	174090
Working pressure	max. 6 bar (stand by 5 bar)	Wall mount	174280
Steam temperature	max. 164°C at 6 bar	Console for splash protection (for wall mounting)	174288
Boiler capacity	2.2 l		
Safety class	IP 33		
Duty	Interrupt operation 5/5 min		

\*Other voltages on request

Technical changes reserved



Wasi-Steam 2 without accessories



Steam gun



Splash protection

### Major features at a glance

- powerful and efficient steam-jet for all cleaning applications
- clear pressure and temperature displays
- simple operation, high safety standards
- durable, versatile and ecofriendly
- stainless steel housing
- water connection, electronic delimer, service interval display and huge stainless steel boiler!

[www.wassermann-dental.com](http://www.wassermann-dental.com)