

● FINE BLASTING UNIT

CEMAT-NT 4/-NT 3/ -NT 2

Surface processing to perfection



Cemat-NT 4

Technical information

Whether sand-blasting, roughening or polishing: first-class work results are guaranteed when you use the precision fine blasting units from the Cemat-NT series. These fine blasting specialists, equipped with anything from 2 to 4 blasting medium containers, combine everything that can be desired of an appliance of this quality: This begins with the even, pulse-free blast and continues from the innovative doser unit via the ultra-rapid start-stop technology right up to the individual and pneumatic blasting control by means of key contact directly on the stylus. The specially coated protection shield provides optimum visibility and especially high protection.

- multifunctional fine blasting appliances with 2 to 4 chambers
- even, pulse-free blast from 1–6 bar
- individual, pneumatic blasting control by means of key contacts directly on the stylus
- practical colour coding system
- high hygiene standard via changeable sleeves
- no electronic components, with the exception of lighting
- filter relief sand drawer for the catchment of coarse materials
- with handy air-expellant jet for simple cleaning

● Technical data and accessories

Cemat-NT 4	Item no.: 190992
Cemat-NT 3	Item no.: 190990*
Cemat-NT 2	Item no.: 190993**
Voltage*	220–240 V / 50/60 Hz
Power consumption	0.2 A
Female connector power consumption max.	8 A
Output	50 W
W x H x D	550 x 480 x 610 mm 430 x 480 x 610 mm ^{*/ **}
Weight	42.6 kg, 37.7 kg*, 34.0 kg ^{**}
Sound level	≤ 70 dB (A)
Tank capacity	800 cm ³ = 1.2 kg
Compressed air connection	3-6 bar
Working pressure	1-6 bar
Extractor nozzle	75 mm

We recommend the use of Wassermann blasting materials and our suction unit SG-10 / SG-20.

Incl. in delivery	Item no.:
PVC hose Ø 8x6 textile (blue), 2 m	320026
Supply line	592012
Cemat-NT 4 including 4 fine blasting jets 0.6 (yellow), 0.8 (blue), 0.8 (red), 1.2 (green)	190992
Cemat-NT 3 including 3 fine blasting jets 0.6 (blue), 0.8 (red), 1.2 (green)	190990
Cemat-NT 2 including 2 fine blasting jets 0.8 (red), 1.2 (green)	190993
Accessories	Item no.:
Automatic stand	190420
PVC hose Ø 8x6 textile (blue), per meter	320026
Sleeve	190175
Spare fine blasting jet 0.4 mm	190254
Spare fine blasting jet 0.6 mm	190250
Spare fine blasting jet 0.8 mm	190251
Spare fine blasting jet 1.2 mm	190252

*Other voltages on request

Technical changes reserved



Cemat-NT 2



Cemat-NT 3



Cemat-NT 4



Major features at a glance

- precision long-life fine blasting units for blasting, roughening and polishing
- permanent pressure system against pollution and dampness in the blasting medium
- individual blasting pressure for each stylus
- problem-free change of blasting medium for all usual blasting media sizes possible over all chambers without use of tools
- even blasting, low consumption
- specially coated protection shield for optimal visibility

www.wassermann-dental.com