

● MODEL TRIMMER HSS-99

The power pack for all occasions



HSS-99 with special accessoires

Technical information

Powerful, innovative and state of the art – the HSS-99 stands for high-end technology and versatile, efficient usage. Powerful 1200 Watt motor suitable for the hardest plasters. Its especially generous disc and maximised working section also allow large models to be trimmed. The top model from the Wassermann plaster trimmer range is the result of years of experience and development based on extensive discussions between the lab and users.

The result: Perfection down to the last detail.

- fully diamond-impregnated grinding disc for silky smooth and flat model surfaces
- high grinding speed of 3000 min⁻¹
- extremely powerful, maintenance-free 1200 Watt motor for high efficiency, splash-protected as per IP 54
- optimal water control using ball valve and solenoid valve
- extra-large work section of 12.5 x 16 cm
- easy to open using rear handle
- powerful disc and interior rinsing

● Technical data and accessories

HSS-99 without grinding disc	Item no.: 164999
Voltage*	220–240 V / 50/60 Hz
Power consumption	5.5 A
Output	1200 W
W x H x D	400 x 430 x 470 mm
Weight	33.0 kg
Sound level	65 dB (A)
Speed	3000 min ⁻¹
Duty	S6-40 %
Table dimensions projecting	160 mm
HSS-99 incl. grinding disc	Artikel Nr.:
incl. diamond grinding disc	164998

Incl. in delivery	Item no.:
Inclined support	164037
Supply hose 1 m	320006
Discharge hose 1 m	320017
Accessories	Item no.:
Spare diamond grinding disc including puller and grease	164063
Grinding device KMS according to B. E. G. G.	164080
Directional light GTE	164081
Template for upper jaw to trim the lower jaw tube level	168082
Sealing grease for trimmer, 25 ml	629002
Puller for grinding disc change	164026

*Other voltages on request

Technical changes reserved



HSS-99 with inclined support



HSS-99 with special accessories for orthodontics



Major features at a glance

- top-class plaster trimmer
- fully diamond-impregnated disc, Ø 300 mm
- 3000 min⁻¹
- 1200 Watt high-performance motor
- numerous expansion possibilities
- high-quality, plastic-powder-coated housing
- adjustable drain supports for easier installation
- expandable by disc equipment and/or adjustable light for 3 D trimming – ideal for orthodontics
- the grinding attachment KMS can be delivered also on swiss- and EOS-version

www.wassermann-dental.com