

“Designed
for superior
results.

simple · powerful · flexible

For curing various materials in dental technology and dentistry to achieve precise results with the materials used. Our DEG-1 pressure compactor compresses the residual air under pressure so effectively that it ensures a bubble-free, homogeneous result.

OUR PROVEN FORMULA FOR SUCCESS:

We combine technology, experience, expertise, and team spirit.



Always here for you

Do you want to request a quote, order products, or have questions about your order, delivery, or warranty?

We look forward to your call and are happy to assist you.

DEG-1
Compactor



Info-Hotline

+49 (0)40 730 926-0



www.wassermann.hamburg



fb.com/WassermannDentalmaschinen

SCAN ME



WASSERMANN
Hamburg · Germany

DEG-1 Compactor

For a bubble-free and
homogeneous material
structure



QUALITY MADE IN GERMANY SINCE 1927

The easy operation and extreme durability reliably provide ease of use. Depending on the inlet pressure, a variable compression pressure of up to 6.0 bar can be achieved, which can then be reduced as needed.



Key Features at a Glance:

- For bubble-free, homogeneous material structure and smooth raw casts
- Integrated stainless steel insert, designed for multiple flasks
- Pressure relief valve for optimal safety
- Non-return valve for consistent operating pressure
- Automatic ventilation: The pressure chamber is automatically vented after the set time
- Included in delivery:
2 m compressed air hose
Power supply

DEG-1
Item No. 180999

Technical Specifications

Operating Pressure: max. 6 bar
Current Consumption: 0.1 A

Dimensions (W x H x D): 270 x 260 x 380 mm
Pressure Chamber (H x D): 130 x 250 mm

Weight: 13 kg
Sound Level: ≤ 70 dB(A)

Info & Consultation:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15 - 17
21031 Hamburg, Germany

Tel: +49 (0)40 730 926-0
www.wassermann-dental.com
sales@wassermann-dental.com