

“developed
for more
flexibility

- efficient & powerful

Our new **Wamix Economy** is a user-friendly vacuum mixer for all types of dental plasters, investments and silicones.

Our proven recipe for success:

We combine technology, experience, competence and team spirit.



Always available

Would you like to request a quote or place an order? Or do you have questions about your order, delivery or warranty?

We look forward to your call and will be happy to help you.

Wamix Economy Vacuum Mixer



Info-Hotline

+49 (0)40 730 926-0



www.wassermann.hamburg



fb.com/WassermannDentalmaschinen

SCAN ME



Wamix Economy Vacuum Mixer

NEU Now simpler and
faster to use


QUALITY MADE IN GERMANY SINCE 1927

The mixing device convinces with good and bubble-free mixing results. The simple operation and extreme robustness reliably make work easier. A newly control system enables the user to work quickly and easily.



The most functions and features at a glance:

- Designed for long life
- Flexible use as wall and table unit (table stand not included)
- Low-maintenance and easy to clean
- For all dental plasters, investments and silicones
- Digital timer
- Mixing bowl and stirrer are included
- Inclined motor for easy attachment of the mixing cup
- Comfortable one-hand operation, optimal mixing results
- High-quality workmanship, made of metal
- Can be mounted in just a few steps

Notes:

Technical Data:

| | | | | | |
|--------------------|--------|-------------------|------------------------|-------------------|--------|
| Output: | 94 W | W x H x D (max.): | 195 x 310 x 180 mm | Weight: | 9,3 kg |
| Power consumption: | 0,44 A | Vacuum: | ca. 120 mbar (absolut) | Protection class: | IP 22 |

Support Service:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15 - 17
21031 Hamburg

Tel: +49 (0)40 730 926-0
www.wassermann-dental.com
sales@wassermann-dental.com