

Recommended Retail Prices 2025



## **PRICE LIST 2025**

This price list is effective from 01.01.2025, ex works Hamburg and is subject to change without prior notice and cannot be considered binding to the seller for any specific period.

Prices are quoted in EUR (EURO), without V.A.T. and subject to our normal terms and conditions of sale regarding delivery and payment. You could request these terms by phone or download them from our homepage.

Specifications are subject to alteration without prior notice.

Orders placed below EUR 50.-- are subject to a handling fee of EUR 10.--.

The plastic coating usually is:  
natural white close to RAL 9016 and black close to RAL 9005.

Colours are also subject to change without notice. Alternative colours are available at additional surcharge.

All units are available in 220-240 V / 50/60 Hz.  
Other voltages upon enquiry. Please let us know.  
Units equipped with heating are only available in 220-240 V / 50/60 Hz.

**Wassermann Dental-Maschinen GmbH**  
**Rudorffweg 15 – 17, 21031 Hamburg, Germany**

Export Department:

Phone: +49 40 730 926 22 or +49 40 730 926 24

[export@wassermann-dental.com](mailto:export@wassermann-dental.com)

[www.wassermann.hamburg](http://www.wassermann.hamburg)

[www.facebook.com/WassermannDentalmaschinen](https://www.facebook.com/WassermannDentalmaschinen)

## Contents

Modell production.....	3
Plaster Dispenser and Accessories.....	3
Dosing scales.....	4
Accessories for DSW-2/ -2/4K.....	4
Vacuum mixers.....	4
Accessories for vacuum mixers.....	5
Vibrators.....	5
Accessories for vibrators.....	6
Trimmers.....	6
Accessories for trimmers.....	7
Spare grinding discs.....	8
Modelling.....	9
Laboratory technology machines.....	9
Accessories for laboratory technology machines.....	10
Steam cleaners.....	10
Accessories Wasi-Steam Classic/ Classic II.....	12
Splash protectors.....	12
Boil-Out and polymerisation units.....	12
Accessories for boil out and polymerisation units.....	13
Compactors.....	14
Mixing and dosing unit.....	14
Accessories for mixing and dosing unit.....	14
Ejection press.....	15
Flask presses.....	15
Accessories for flask presses.....	16
Surface finish.....	16
Fine blasting units.....	16
Accessories for Gobi-2 (Item no. 190896).....	19
Accessories for Puna (Item no. 190897).....	19
Polishing compact units.....	19
Polishing motor.....	20
Polishing troughs.....	20
Accessories for polishing lathes.....	20
Filter for polishing compact units.....	21
Automatic stand.....	21
Suction technology.....	21
Suction Box.....	23
Accessories for suction units.....	23
Filters for suction units.....	24
Consumable materials.....	25
Duplicating silicone.....	25
Blasting materials.....	26

Price in €

**Modell production****Plaster Dispenser and Accessories**

<b>166997 D</b>	<b>Plaster Dispenser WGS-25W II</b> Wall Model, rotatable (180 °), for 25 kg plaster container	4.575,00
<b>166115 D</b>	<b>Table-Top Support</b> for WGS-25 W II	658,00

**Dosing scales**

<b>172988 A</b>	<b>DSW-2 (with 1 pump, not cooled)</b> For exact dosage of plaster and investment material You can order up to 5 additional pumps	2.460,00
-----------------	---	----------

**Accessories for DSW-2/ -2/4K**

<b>172299 C</b>	<b>Additional pump equipment (not cooled)</b>	224,00
-----------------	---	--------

	<b>Price in €</b>
<b><u>Vacuum mixers</u></b>	
<b>171970 A Wamix-Touch</b> with vacuum pump, incl. wall holder and mixing cup 350 ml with touch screen for fast and easy operation	2.205,00
<b>171971 A Wamix-Touch Injector</b> with injector technique, incl. wall holder and mixing cup 350 ml with touch screen for fast and easy operation	1.990,00
<b>171972 B Wamix-Economy</b> with vacuum pump, incl. wall holder and mixing cup 350 ml	1.435,00

**Accessories for vacuum mixers**

**for Wamix-Touch/ Wamix-Touch Injector and  
Wamix-2/ -2M/ -2MS/ -3M/ -2000/ CCI<sup>1</sup>**

**CCI<sup>1</sup>** = Wamix-Comfort Plus, -Comfort, -Classic, -Injector

<b>171295 C</b> Mixing cup Wamix mini, 75 ml (75 ml mixing cup not suitable for beaker mount Wamix-Comfort Plus)	215,00
<b>171230 C</b> Mixing cup Wamix small, 200 ml	224,00
<b>171283 C</b> Mixing cup Wamix medium, 350 ml	231,00
<b>171237 C</b> Mixing cup Wamix large, 600 ml	239,00
<b>171269 C</b> Mixing cup Wamix extra large, 1000 ml (1000 ml mixing cup not suitable for beaker mount Wamix-Comfort Plus)	305,00
<b>171744 C</b> Table stand Wamix-Touch/ Wamix-Touch Injector	350,00
<b>171122 C</b> Set of sealings (set = 5 pcs.)	8,80
<b>610025 C</b> Sinterfilter	24,00

	<b>Price in €</b>
<b><u>Vibrators</u></b>	
<b>169997 A</b> <b>KV-56</b> , <u>plastic-powder-coated</u> cast housing White Table dimensions 430 x 286 mm, vibration intensity infinitely variable	1.650,00
<b>169999 A</b> <b>KV-36</b> , <u>plastic-powder-coated</u> stainless steel White Table dimensions 293 x 200 mm, vibration intensity infinitely variable	795,00
<b>169992 A</b> <b>KV-26 plus</b> , <u>plastic-powder-coated</u> stainless steel White Table dimensions 210 x 130 mm, vibration intensity infinitely variable	689,00
<b>169996 A</b> <b>KV-26</b> , <u>plastic-powder-coated</u> stainless steel White Table dimensions 210 x 130 mm	565,00
<b>169998 A</b> <b>KV-16</b> , <u>plastic-powder-coated</u> stainless steel White Table dimensions 130 x 130 mm	490,00
<b>169989 A</b> <b>KV-36</b> , <u>matt</u> stainless steel Table dimensions 293 x 200 mm, vibration intensity infinitely variable	795,00
<b>169988 A</b> <b>KV-26 plus</b> , <u>matt</u> stainless steel Table dimensions 210 x 130 mm, vibration intensity infinitely variable	689,00
<b>169987 A</b> <b>KV-26</b> , <u>matt</u> stainless steel Table dimensions 210 x 130 mm	565,00
<b>169986 A</b> <b>KV-16</b> , <u>matt</u> stainless steel Table dimensions 130 x 130 mm	490,00
<b><u>Accessories for vibrators</u></b>	
<b>169510 C</b> Spare rubber plate KV-56	136,00
<b>169010 C</b> Spare rubber plate KV-36	64,00
<b>169101 C</b> Spare rubber plate KV-26/ KV-26 plus	55,50
<b>169201 C</b> Spare rubber plate KV-16	31,50

	<b>Price in €</b>
<b><u>Trimmers</u></b>	
<b>164998 A HSS-99 - 1200 Watt – incl. diamond grinding disc</b> grinding disc entirely diamonded, directly coated, incl. solenoid valve and ball tap, large grinding area	4.885,00
<b>168988 A HSS-88 Super 1.500 Watt incl. Diamond grinding disc</b> incl. solenoid valve, grinding disc entirely diamonded, directly coated	3.990,00
<b>168998 A HSS-88 - 1000 Watt - incl. diamond grinding disc</b> incl. solenoid valve, grinding disc entirely diamonded, directly coated	3.845,00
<b>168996 A HSS-88 - 1000 Watt - incl. Dia-Quick grinding disc</b> incl. equipment for first installation and solenoid valve, grinding disc partly diamonded, directly coated	3.385,00
<b>163990 B HSS-TR Dry Trimmer Dia-Quick</b> with partially diamond disc	3.235,00
<b>168989 A HSS-88 Super 1.500 Watt incl. corundum-grinding disc</b> incl. solenoid valve, corundum-grinding disc 10 mm thick	2.535,00
<b>168999 A HSS-88 - 1000 Watt – incl. corundum-grinding disc</b> incl. solenoid valve, corundum-grinding disc 10 mm thick	2.345,00
<b>168991 A HSS-88 - 1000 Watt - incl. foil grinding disc</b> incl. solenoid valve, incl. equipment for first installation	2.345,00
<b>162997 A HSS-AZ - 750 Watt – incl. Dia-Quick grinding disc</b> incl. solenoid valve, grinding disc partly diamonded, directly coated	2.615,00
<b>162996 A HSS-AZ - 750 Watt - incl. foil grinding disc</b> incl. solenoid valve	1.640,00
<b>161996 A HSS-ZA</b> <b>incl. Dia-Quick grinding disc and brush spindle on left side</b> combination trimmer/polishing motor -500 Watt-	3.665,00
<b>161997 A HSS-ZA</b> <b>incl. foil grinding disc and brush spindle on left side</b> combination trimmer/polishing motor -500 Watt-	2.555,00
<b>161114 C Polishing trough base for HSS-ZA (2017)</b>	80,50

(Polishing troughs see page 19)

**Price in €**

**Accessories for trimmers**

**Grinding attachment for HSS-99**

**for producing 3D orthodontic models:**

<b>164130 D</b>	Model table KMS according to Begg	2.230,00
<b>164110 D</b>	Laser light guide for model table	870,00
<b>164034 D</b>	Orthodontic template, lower jaw/ upper jaw/ tuber plane	54,00

**Grinding attachment for HSS-88**

**for producing 3D orthodontic models:**

<b>168140 D</b>	Model table KMS according to Begg	2.230,00
<b>168120 D</b>	Laser light guide for model table	911,00
<b>164034 D</b>	Orthodontic template, lower jaw/ upper jaw/ tuber plane	54,00

**!** The Grinding attachment KMS can be delivered also as Swiss and EOS-versions. Price on demand. **!**

<b>164161 C</b>	Tray for trimmer HSS-99 fitting in a sink (We always need dimensions of the sink and design 1:1)	329,00
<b>168160 C</b>	Tray for trimmer HSS-88/ AZ fitting in a sink (We always need dimension of the sink and design 1:1)	284,00
<b>164037 C</b>	Inclined bearing table for HSS-99	105,00
<b>168037 C</b>	Inclined bearing table for HSS-88/HSS-AZ/HSS-ZA	103,00
<b>168168 C</b>	Model tray for HSS-88	127,00
<b>168166 C</b>	Higher base for trimmers	147,00
<b>320006 C</b>	1 m water supply hose ½", with cloth inlet and 2 hose clips	25,00
<b>320017 C</b>	1 m drainage hose 1", translucent, made of plastic	31,50
<b>629002 C</b>	Sealing lubricant for trimmer, 25 ml	8,80
<b>164026 C</b>	Puller for grinding disc exchange for HSS-99/ HSS-88/ HSS-AZ/ HSS-ZA	45,50



	<b>Price in €</b>
<b><u>Spare grinding discs</u></b>	
<b>164063 C Spare diamond grinding disc for HSS-99,</b> entirely diamonded, directly coated, Ø 300 mm, incl. puller and sealing lubricant	2.330,00
<b>168063 C Spare diamond grinding disc for HSS-88,</b> entirely diamonded, directly coated, Ø 260 mm, incl. puller and sealing lubricant	1.980,00
<b>162063 C Spare Dia-Quick grinding disc for HSS-AZ/-88/-ZA,</b> partly diamonded, directly coated, Ø 255 mm, incl. sealing lubricant	1.145,00
<b>168024 C Spare corundum grinding disc for HSS-88,</b> 10 mm thick, Ø 260 mm, incl. puller and sealing lubricant	356,00
<b>160020 C Spare corundum grinding disc for HSS-65/HSS-68,</b> 10 mm thick, Ø 242 mm, incl. puller and sealing lubricant	390,00
<b>162013 C Spare foil grinding disc for HSS-AZ/-88/-ZA,</b> 1 set = 4 pieces, Ø 250 mm, incl. sealing lubricant	178,00
<b>162023 C Spare foil grinding disc (fine) for HSS-AZ/-88/-ZA,</b> 1 set = 4 pieces, Ø 250 mm, incl. sealing lubricant	216,00
<b>168071 C Conversion kit for HSS 88</b> for Dia-Quick and foil grinding disc incl. puller for grinding disc change and sealing lubricant	336,00

Price in €

**Modelling****Laboratory technology machines****HSM-50 / W&H Perfecta 600**

1000 – 50000 min<sup>-1</sup>, 200 Watt, infinitely adjustable motor,  
 engine torque 7,8 N/cm, cpl. with rest for handpiece

<b>111986 A</b>	<b>Motor handpiece with foot-operated starter</b> with blow-out function operating controls as accessories available	2.460,00
<b>111987 A</b>	<b>Motor handpiece with knee-operated starter</b> with blow-out function incl. operating controls foot-operated starter with blow-out function as accessories available	2.720,00
<b>111988 A</b>	<b>Motor handpiece with table-operated starter</b> incl. operating controls and foot-operated starter with blow-out function	2.720,00

**HSM-40 / W&H Perfecta 300**

1000 – 40000 min<sup>-1</sup>, 160 Watt, infinitely adjustable motor,  
 engine torque 5,5 N/cm, cpl. with rest for handpiece

<b>111980 A</b>	<b>Motor handpiece with foot-operated starter</b> operating controls as accessories available	1.830,00
<b>111981 A</b>	<b>Motor handpiece with knee-operated starter</b> incl. operating controls foot-operated starter as accessories available	2.125,00
<b>111982 A</b>	<b>Motor handpiece with table-operated starter</b> incl. operating controls and foot-operated starter	2.125,00

**Accessories for laboratory technology machines**

<b>111121 C</b>	<b>Motor handpiece 50 / 40 / 300 / 600 (without supply line)</b>	1.590,00
<b>111117 C</b>	<b>Operating controls HSM-50 / W&amp;H Perfecta 600</b>	344,00
<b>111111 C</b>	<b>Operating controls HSM-40 / W&amp;H Perfecta 300</b>	344,00
<b>111124 C</b>	<b>Holder for operating controls 50 / 40 / 300 / 600</b>	49,50
<b>111119 C</b>	<b>Foot starter for 50* / 40 / 300 / 600*</b> for knee-/ and table-operated starter	229,00

\*HSM-50 / W&amp;H Perfecta 600 with blow out function

Price in €

Steam cleaners
**Wasi-Steam Compact and Wasi-Steam 3  
 (fixed water connection required)**

Thorough and efficient steam cleaner suitable for all cleaning jobs in the lab. Choice of wet or dry steam

**174989 B Wasi-Steam 3 Stainless Steel New for the IDS!** 3.590,00

174994 A *Wasi-Steam Compact Stainless Steel White, remaining stock!* 3.590,00

174994 A *Wasi-Steam Compact Stainless Steel, remaining stock!* 3.590,00

174987 A *Wasi-Steam Compact Stainless Steel with filling container* 3.950,00

**174995 C Wasi-Clean Water softener (strongly recommended!)** 292,00

Table top unit W 190 x H 405 x D 260 mm

**Wasi-Steam Classic II (manual filling)**

Compact steam cleaner with steam gun, manual filling

Flexible choice of location. Operate with distilled/softened water recommended!

**174996 A #CCWH# Wasi-Steam Classic II Stainless Steel White** 1.695,00

**174996 A #CCVA# Wasi-Steam Classic II Stainless Steel** 1.695,00

174970 C *Wasi-Steam Box #CCWH# Stainl. Steel White, remaining stock!* 2.040,00  
*Compact working cabin with integrated suction*

**174970 C Wasi-Steam Box #CCVA# Stainless Steel** 2.040,00  
 Compact working cabin with integrated suction

	<b>Price in €</b>
<b>174633 C Quantity Control Set for Wasi-Steam 3</b> To regulate <u>wet steam saturation</u> and <u>steam quantity</u> . It would require to retrofit the unit at our factory or at one of our authorized service stations.	230,00
<b>174414 C Quantity Control Set for Wasi-Steam Compact</b> To regulate <u>wet steam saturation</u> and <u>steam quantity</u> . It would require to retrofit the unit at our factory or at one of our authorized service stations.	230,00
 <b><u>Accessories Wasi-Steam Classic/ Classic II</u></b> 	
<b>655002 C Lime remover, 1000 ml</b>	103,00
<b>174358 C Set for wall mount</b>	22,50
<b>390100 C Tank Closing Cap WSC</b>	29,50
 <b><u>Splash protectors</u></b> 	
<b>174460 C Splash protection II</b> W 250 x H 265 x D 255 mm open at the front and top entirely made of stainless steel	245,00
<b>174360 C Splash protection</b> W 235 x H 385 x D 294 mm top front with fixed stainless steel hood (hood: open top, chimney effect) with plastic basement at the bottom	245,00
<b>174230 C Splash protection stainless steel Wasi-Steam 3</b>	389,00

	<b>Price in €</b>
<b><u>Boil-Out and polymerisation units</u></b>	
<b>170970 D Wapo-Ex 12 II</b> Boil out unit for 12 flask-halves, complete with baskets, equipped with touch screen for fast and easy operation, incl. lime remover and wax remover	9.215,00
<b>170977 C Wapo-Ex 8 Plus</b> (combination automat) Compact tabletop automat for boil out up to 8 flask halves or polymerisation up to 8 flasks in a clamp, equipped with touch screen for fast and easy operation, complete with baskets and incl. service set	5.125,00
<b>170980 C Wapo-Mat III</b> Polymerisation unit for longtime- and boil polymerisation, for up to 9 flasks in the strainer, equipped with touch screen for fast and easy operation, incl. lime remover	2.360,00
<b>170993 C Polyquick 260-T105</b> (combination automat) Pressure polymerisation unit for hot and cold polymerisation, pressure polymerisation up to 105 °C in a water bath, can also be used as a pressure compactor, very large capacity Ø 240 mm, usable depth 160 mm, suitable for extra large articulators, touch screen for fast and easy operation	4.470,00
<b>170976 C Polyquick 160-T95</b> (combination automat) Pressure polymerisation unit for hot and cold polymerisation, pressure polymerisation up to 95 °C in a water bath can also be used as a pressure compactor, capacity Ø 150 mm, depth 285 mm touch screen for fast and easy operation Optional feature: Dry pressure compression up to approx. 70°C	3.590,00

	<b>Price in €</b>
<b><u>Accessories for boil out and polymerisation units</u></b>	
<b>170120 C</b> Moveable stand for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 with joint rollers, W x H x D: 620 x 550 x 430 mm	685,00
<b>170100 C</b> Moveable stand for Wapo-Ex with joint rollers, W x H x D: 560 x 550 x 410 mm	685,00
<b>170330 C</b> Service set for Wapo-Ex 12 II/ Wapo-Ex 12 consisting of lime remover and wax remover	218,00
<b>170021 C</b> Service set for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 consisting of 4 filter mats, lime remover, wax remover	267,00
<b>170020 C</b> Service set for Wapo-Ex consisting of 10 filter mats, lime remover, wax remover	267,00
<b>611007 C</b> Spare filter for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 / Wapo-Ex old up to #40189 consisting of 4 filter mats 400 x 325 mm	64,00
<b>611013 C</b> Spare filter for Wapo-Ex consisting of 10 filter mats 400 x 145 mm	64,00
<b>655001 C</b> Wax remover, 2 liter	127,00
<b>655002 C</b> Lime remover, 1000 ml	103,00
<b>!</b> When you order new filters, please mention the S/N serial number of the unit.	<b>!</b>

### **Compactors**

<b>170993 C</b> Polyquick 260-T105 (combination automat, see page 12)	4.470,00
<b>170976 C</b> Polyquick 160-T95 (combination automat, see page 12)	3.590,00
<b>180994 C</b> <b>Compactor DEG-1 New Design!</b> old item no 180999	2.450,00
<b>380032 C</b> O-Ring for DEG / Polyquick	31,35
<b>380063 C</b> O-Ring for Polyquick 260	82,00

Price in €

Mixing and dosing unit

<b>172970 C</b>	<b>Sidomix III New Design!</b>	1.950,00
	Dosage and mixing unit for addition-vulcanizing duplicating silicone, 1:1	

Accessories for mixing and dosing unit

<b>172640 C</b>	<b>Set material cask, cpl. 2 x 2 liter incl. quick fitting coupling</b>	116,00
<b>172633 C</b>	<b>Connecting tubes cpl., 2 pcs.</b>	75,50
<b>172634 C</b>	<b>Adapter set for WaSil silicone 2 x 1 kg cask</b>	68,00
<b>172636 C</b>	<b>Adapter set for WaSil silicone 2 x 6 kg cask</b>	75,50
<b>622002 C</b>	<b>Mixing canulas 1:1, 5.4 mm inner Ø, pink, 100 pcs. connectable to dosage automat 1:1 Sidomix II, Sidomix, DA 2000 and similar equipment</b>	85,50

**! For a long durability of the mixing and dosing unit we recommend the use of our Wassermann WaSil duplicating silicone. More on Page 26. !**

	<b>Price in €</b>
<b><u>Ejection press</u></b>	
<b>128999 A APR-1</b> for cast rings/ cast flasks (DeguDent/ Heraeus Kulzer), cast ring set 1	1.155,00
<b><u>Flask presses</u></b>	
<b>126999 A Hydraulic press WW-33 for 1-3 flasks</b>	1.490,00
<b>126998 A Hydraulic press WW-33 for 2-4 flasks</b>	1.610,00
<b>129997 C Universal flask clamp for 1 flask</b>	217,00
<b>129996 C Universal flask clamp for 2 flasks</b>	224,00
<b>129995 C Universal flask clamp for 3 flasks</b>	229,00
<b><u>Accessories for flask presses</u></b>	
<b>129994 C Stainless steel clamp clip for 1 flask</b>	163,00
<b>129993 C Stainless steel clamp clip for 2 flasks</b>	170,00
<b>129992 C Stainless steel clamp clip for 3 flasks</b>	180,00
<b>129990 C Universal flask</b>	141,00
<b>629005 C Spare bottle hydraulic oil, 500 ml</b> for WW-11, WW-22, WW-33 and KÜVOMAT	103,00
<b>629004 C Hydraulic Oil, Bottle 100 ml</b>	36,50



Price in €

**Surface finish****Fine blasting units**

<b>190986 A Cemat-4 Plus with intuitive hand controls</b>	7.515,00
with 4 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	
<b>190984 A Cemat-4 with 4 blasting material containers and sand-drawer, individually adjustable working pressure per chamber</b>	6.450,00
<b>190985 A Cemat-3 Plus with intuitive hand controls</b>	6.315,00
with 3 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	
<b>190983 A Cemat-3</b>	5.795,00
with 3 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	
<b>190981 A Cemat-2 Plus with intuitive hand controls</b>	5.675,00
with 2 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	
<b>190982 A Cemat-2</b>	5.155,00
with 2 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	
<b>190884 A Gobi Comfort-4</b>	4.820,00
with 4 blasting material containers and sand-drawer,	
<b>190883 A Gobi Comfort-3</b>	4.205,00
with 3 blasting material containers and sand-drawer,	
<b>190882 A Gobi Comfort-2</b>	3.540,00
with 2 blasting material containers and sand-drawer,	
<b>190881 A Gobi Classic</b>	2.460,00
with 2 blasting material containers	
<b>190885 A Puna II</b>	1.695,00
with 2 blasting material containers	

! The Wassermann fine blasting units are optimally designed for use with our Wassermann dental extraction systems. Page 22. !

For optimal durability, we recommend Wassermann blasting material to extend the service life of our units. Page 27.

Price in €

**Accessories for fine blasting units**
**Accessories Cemat II/NT, Gobi Comfort/Classic, Puna/II**

For Cemat OLD please see page 20

<b>190553 C</b>	Spare fine blasting nozzle 0.4 mm (only Cemat-NT)	110,00
<b>190550 C</b>	Spare fine blasting nozzle 0.6 mm	110,00
<b>190551 C</b>	Spare fine blasting nozzle 0.8 mm	110,00
<b>190552 C</b>	Spare fine blasting nozzle 1.2 mm	110,00
<b>190565 C</b>	Spare nozzle (tank) 2.0 mm	18,50
<b>190566 C</b>	Spare nozzle (tank) 1.4 mm	18,50
<b>190494 C</b>	Key for assembly/ disassembly of the tank nozzle	33,50
<b>190533 C</b>	Rep. Set Sleeve Gobi Cemat II	70,50

**Accessories for Cemat II**

<b>190489 C</b>	Spare hose -yellow-	109,00
<b>190485 C</b>	Spare hose -blue-	109,00
<b>190475 C</b>	Spare hose -red-	109,00
<b>190487 C</b>	Spare hose -green-	109,00
<b>190432 C</b>	Spare protection sheet with special coating	110,00

**Accessories for Gobi Comfort and Gobi Classic**

<b>190739 C</b>	Spare hose -yellow-	103,00
<b>190738 C</b>	Spare hose -blue-	103,00
<b>190729 C</b>	Spare hose -red-	103,00
<b>190728 C</b>	Spare hose -green-	103,00
<b>190432 C</b>	Spare protection sheet with special coating	110,00

**Accessories for Puna II**

<b>190475 C</b>	Spare hose -red-	109,00
<b>190485 C</b>	Spare hose -blue-	109,00
<b>190686 C</b>	Spare protection sheet with special coating	110,00

Price in €

**Accessories for fine blasting units (previous versions)**
**Accessories for Cemat-NT  
 (Item no. 190992/ 190990/ 190993)**

Spare Nozzles see page 18

<b>190432 C</b>	Spare protection sheet with special coating from serial number 128538	110,00
<b>190458 C</b>	Conversion kit protection sheet with special coating up to serial number 128537	126,00
<b>190230 C</b>	Spare protection sheet up to serial number 128537	110,00
<b>612003 C</b>	Spare protection foil (plastic), 1 set = 5 pieces, up to serial number 128537	83,50
<b>190239 C</b>	Spare hoses	126,00
<b>190533 C</b>	Rep. Set Sleeve Gobi Cemat II	70,50

**Accessories for Gobi-2 (Item no. 190896)**

Spare Nozzles see page 18

<b>190432 C</b>	Spare protection sheet with special coating from serial number 128423	110,00
<b>190636 C</b>	Conversion-kit protection sheet with special coating # 124048 - # 127319	397,00
<b>612004 C</b>	Spare protection foil (plastic), 1 set = 5 pieces, up to serial number 127319	83,50
<b>190487 C</b>	Spare hose for green chamber	109,00
<b>190475 C</b>	Spare hose for red chamber	109,00
<b>190485 C</b>	Spare hose for blue chamber	109,00
<b>190489 C</b>	Spare hose for yellow chamber	109,00
<b>190533 C</b>	Rep. Set Sleeve Gobi Cemat II	70,50

**Accessories for Puna (Item no. 190897)**

Spare Nozzles see page 18

<b>190686 C</b>	Spare protection sheet with special coating	110,00
<b>190475 C</b>	Spare hose for red chamber	109,00
<b>190485 C</b>	Spare hose for blue chamber	109,00
<b>190533 C</b>	Rep. Set Sleeve Gobi Cemat II	70,50

**Accessories for Cemat OLD  
 (Item no. 190996/ 190999/ 190998/ 190997)**

<b>190080 C</b>	Spare fine blasting nozzle 0.4 mm	126,00
<b>190081 C</b>	Spare fine blasting nozzle 0.6 mm	126,00
<b>190083 C</b>	Spare fine blasting nozzle 0.8 mm	126,00
<b>190084 C</b>	Spare fine blasting nozzle 1.2 mm	126,00
<b>190111 C</b>	Spare hoses	103,00
<b>612001 C</b>	Spare protection foil, 1 set = 5 pieces	83,50

	<b>Price in €</b>
<b><u>Polishing compact units</u></b>	
<b>146979 B Polishing unit WP-Ex 3000 II</b>	5.740,00
incl. adjustable motor (750-6000 min <sup>-1</sup> ), dust extraction, LED illumination and LED Spot, rubber polishing troughs, drawer, brush spindle left + right, filters, 2 storage bowls	
<b>146978 B Polishing unit WP-Ex 2000 II</b>	4.675,00
incl. 500 W motor, dust extraction, LED illumination, rubber polishing troughs, drawer, brush spindle left + right, filters, 2 storage bowls	
<b>146971 B Polishing unit WP-Ex 10 II</b>	3.065,00
incl. 500 W motor, dust extraction, LED illumination, rubber polishing trough, brush spindle right, filter, 1 storage bowl	
<b><u>Polishing motor</u></b>	
<b>155999 B WSM-2 (1500/ 3000 min<sup>-1</sup>)</b>	1.105,00
incl. brush spindle left + right	
<b><u>Polishing troughs</u></b>	
<b>143123 C Polishing trough without extractor nozzle (2017)</b>	266,00
incl. adjustable protection sheet and removable rubber trough	
<b>143115 C Polishing trough with extractor nozzle (2017)</b>	299,00
(Ø75 mm), incl. cap incl. adjustable protection sheet and removable rubber trough	
<b>143131 C Polishing trough with extractor nozzle (2017)</b>	299,00
(Ø40 mm), incl. cap incl. adjustable protection sheet and removable rubber trough	
<b>143550 C Extraction traverse II, with switch to guide the air</b>	597,00
output nozzle Ø75 mm	

Price in €

**Accessories for polishing lathes**

<b>152025 C</b>	Brush spindle 92 mm for left side	55,50
<b>152022 C</b>	Brush spindle 92 mm for right side	55,50
<b>152057 C</b>	Stone holder for left side	55,50
<b>152056 C</b>	Stone holder for right side	55,50
<b>152044 C</b>	Brush spindle 125 mm for left side	83,50
<b>152043 C</b>	Brush spindle 125 mm for right side	83,50
<b>152049 C</b>	Chuck clamping from 1-6 mm	218,00
<b>155123 C</b>	Wrench for pull nut	67,50
<b>146390 C</b>	Daylight LED Spot	301,00
<b>146585 C</b>	<b>RepairSet Glass Panel from old to new</b>	470,00

**Filter for polishing compact units**

<b>611005 C</b>	<b>Spare filter bag for WP-Ex 3000 II/ 2000 II/ 10 II</b> dust class M, for micro fine dust, 1 set = 2 pieces	98,00
<b>611011 C</b>	<b>Spare filter bag WP-Ex 3000/ 2000/ 10</b> dust class M, for micro fine dust, 1 set = 2 pieces	83,50
<b>611021 C</b>	<b>Spare filter bag for WP-Ex 80</b> dust class M, for micro fine dust, 1 set = 2 pieces	83,50

**! Please sent us the serial number of your unit in case of ordering filter material. The serial number allows us to find quickly without any doubt the correct sparteparts of your unit. The identification plate is located on the back or right side of your unit. !**

**Automatic stand**

<b>146170 D</b>	<b>Automatic stand for WP-Ex 3000 II/ 2000 II</b>	3.065,00
	vertically adjustable from 74 – 109 cm Dimension (table top): D 560 x W 675 mm	

	<b>Price in €</b>
<b><u>Suction technology</u></b>	
<b>140997 A Single-place-dust extraction SG-1/2 D (Suction port 75 mm)</b> Brushless motor. Equipped with automatic switch-off when the filter is full. Filter included. Without accessories.	1.885,00
<b>140995 A Single place dust extraction SG-1/1 (Suction port 75 mm)</b> Brushless motor. Filter included. Without accessories.	1.130,00
<b>140991 C Single place dust extraction SG-1/1 Comfort (Suction port 75 mm)</b> Brushless motor. With automatic fill level monitoring. Equipped with auto-start-system.	1.425,00
<b>140970 C Suction unit SG-Economy (Suction port 40 mm)</b> incl. 3 stage filter system, without accessories. <b>Recommended for HSS-TR (dry trimmer)</b> Powerful brushed motor for reliable performance. Equipped with auto-start-system.	1.385,00

	<b>Price in €</b>
<b>145997 A Dust extraction drawer LSG-02 DE</b> complete with universal mounting frame and automatic switch-off in case of exhausted filter capacity, incl. filter , additional equipped with auto-start-system	2.390,00
<b>145998 A Dust extraction drawer LSG-02 D</b> complete with universal mounting frame and automatic switch-off in case of exhausted filter capacity, incl. filter	2.205,00
<b>145240 C Activated charcoal filters for monomer for LSG-02 DE/ LSG-02 D</b> for elimination of dangerous vapors	702,00
<b>145200 C Universal extractor hood for LSG-02 D and DE</b> with exchangeable protection sheet	420,00
<b>145005 C Additional drawer for LSG-02 D and DE</b> can also be retrofitted	227,00
<b>149999 A Mobile dust extraction Timobil</b> cpl. with handpiece motor 1.000 – 40.000 min <sup>-1</sup> , 160 Watt, infinitely adjustable, standard height 900 mm	4.675,00
<b>149130 C Extractor hood for Timobil</b> Closed, can also be retrofitted	401,00
<b>149140 C Extractor hood illuminated for Timobil</b>	629,00
<b><u>Suction Box</u></b>	
<b>143390 C Suction Box Compact III</b> with daylight LED illumination and 1 extractor nozzle Ø 75 mm exterior (optional mounting in the front or in the rear)	805,00
<b><u>Accessories for Suction Box Compact III + previous versions</u></b>	
<b>143303 C</b> Set of arm rests	85,50
<b>143305 C</b> Set of cover plates	221,00
<b>190686 C</b> Spare protection sheet with special coating	110,00

Price in €

**Accessories for suction units****Telescopic channels + corresponding suction parts**

<b>143070 C</b>	Telescopic channel 570 mm	152,00
<b>143090 C</b>	Telescopic channel 455 mm	162,00
<b>143074 C</b>	Telescopic channel 480 mm	161,00
<b>143066 C</b>	Suction part, transverse	120,00
<b>143080 C</b>	Suction part, straight	132,00
<b>145245 C</b>	Activated charcoal filters for monomer - suction part, straight for elimination of dangerous vapors	702,00
<b>143108 C</b>	Cap for suction part transverse	50,00
<b>143200 C</b>	Universal hood for suction part straight	420,00
<b>143220 C</b>	Suction part with wooden rest	350,00
<b>143221 C</b>	Suction part with wooden rest/plastic screen	422,00
<b>143231 C</b>	Suction part with wooden rest/glass screen	470,00

**Special channel + appropriate suction parts**

(e.g. for installation in Freuding tables)

<b>143093 C</b>	Special channel with hose coupling	183,00
<b>143041 C</b>	Extractor nozzle	221,00

**! Channel and one suction part are always both required!**

**!****Hoses**

<b>140057 C</b>	PVC-hose 75 Ø - 1.20 m length, incl. 2 hose clips	110,00
<b>320011 C</b>	PVC-hose 75 Ø - 1.80 m length	156,00
<b>320010 C</b>	PVC-hose 75 Ø - 2.40 m length	192,00
<b>320012 C</b>	PVC-hose 75 Ø - 3.60 m length	269,00
<b>322003 C</b>	Hose clip 75 Ø	5,30
<b>327023 C</b>	PVC-hose 51 Ø inner - per meter	60,00

**Special accessories**

<b>140056 C</b>	Throttle valve 75 Ø	167,00
-----------------	---------------------	--------

**! Reduction- and transition pieces of all measurements available.**

**!**

**Would you please be so kind as to mention inner- and outer diameter in case of supplementary order.**



Price in €

**Filters for suction units****Filters & accessories for suction unit SG-10**

<b>611010 C</b>	Spare filter bag yellow dust class M, for micro fine dust, 1 set = 4 pieces	196,00
<b>611018 C</b>	Spare filter bag white, for plaster class 1 and 2 only (!) regenerative, 1 set = 4 pieces	196,00
<b>610073 C</b>	*Filter bag white (plaster class 1/ class 2), big	77,50
<b>610074 C</b>	*Filter bag yellow (micro fine dust), big	77,50
	<i>* to be used only in combination with item no. 140410</i>	
<b>140410 C</b>	Filter extension module SG-10	271,00
<b>140435 C</b>	Filter cassette, cpl. SG-10/ SG-20	91,00
<b>140424 C</b>	Module activated carbon filter SG-10	539,00
<b>140436 C</b>	Activated carbon filter (AKM 6-ply) SG-10	154,00

**SG-1/1 + SG-2/2 + Timobil**

<b>611010 C</b>	Spare filter bag yellow, dust class M, for micro fine dust, 1 set = 4 pieces	196,00
<b>611018 C</b>	Spare filter bag white, (class 1/class 2) for plaster, regenerable, 1 set = 4 pieces	196,00

**SG-1/2 D+ LSG-02 DE/ LSG-02 D**

<b>611020 C</b>	Spare filter bag yellow, dust class M, for micro fine dust, 1 set = 4 pieces	235,00
<b>611017 C</b>	Spare filter bag white, (class 1/class 2) for plaster, regenerable, 1 set = 3 pieces	178,00
<b>610056 C</b>	Spare filter monomer for LSG-02	71,00

<b>611022 C</b>	<b><u>Filter-Kit SG-Economy</u></b> 1 x Spare-main filter / cartridge filter SG-Economy 5 x Spare-filter bag SG-Economy 2 x Spare-HEPA exhaust filter SG-Economy	231,00
-----------------	---	--------

**Filters & accessories for suction unit SG-20**

<b>610073 C</b>	Filter bag white (plaster class 1/ class 2)	77,50
<b>610074 C</b>	Filter bag yellow (micro fine dust)	77,50
<b>140551 C</b>	Activated carbon filter, cpl. SG-20	149,00
<b>140435 C</b>	Filter cassette, cpl. SG-10/ SG-20	91,00
<b>610019 C</b>	Motor protective filter SG-20	12,00

**WTS-1R**

<b>611019 C</b>	filter bag FB 25, 1 Set = 3 pieces	110,00
-----------------	------------------------------------	--------

Price in €

**Consumable materials****Duplicating silicone****WaSil soft**

Addition-vulcanizing duplicating silicone, 1:1  
 hardness Shore A: 16-18  
 component A white (catalyst)  
 component B yellow (base)

<b>656060 C</b>	<b>WaSil soft 2 x 1 kg component A/ B</b>	116,00
<b>656061 C</b>	<b>WaSil soft 2 x 6 kg component A/ B</b>	594,00
<b>656062 C</b>	<b>WaSil soft 2 x 25 kg component A/ B</b>	2.350,00

**WaSil speed**

Addition-vulcanizing duplicating silicone, 1:1  
 high density, hardness Shore A: > 24  
 component A white (catalyst)  
 component B turquoise (base)

<b>656070 C</b>	<b>WaSil speed 2 x 1 kg component A/ B</b>	132,00
<b>656071 C</b>	<b>WaSil speed 2 x 6 kg component A/ B</b>	731,00
<b>656072 C</b>	<b>WaSil speed 2 x 25 kg component A/ B</b>	2.780,00

**Mixing canulas, 1:1  
 5.4 mm inner Ø**

<b>622002 C</b>	<b>mixing canulas, pink, 100 pcs.</b> connectable to dosage automat 1:1 Sidomix II, Sidomix, DA 2000 and similar equipment	85,50
-----------------	--	-------

Price in €

**Blasting materials****Aluminium oxyde**

For blasting of all dental alloys. The purity, sharp-edged grains and extreme hardness guarantee very efficient blasting.

Silicosis-free, made in Germany

<b>656010 C</b>	<b>Aluminium oxyde 50 µm – 10 kg</b>	95,50
<b>656015 C</b>	<b>Aluminium oxyde 50 µm – 20 kg</b>	197,00
<b>656013 C</b>	<b>Aluminium oxyde 110 µm - 10 kg</b>	95,50
<b>656016 C</b>	<b>Aluminium oxyde 110 µm - 20 kg</b>	197,00
<b>656014 C</b>	<b>Aluminium oxyde 250 µm - 10 kg</b>	95,50
<b>656018 C</b>	<b>Aluminium oxyde 250 µm – 20 kg</b>	197,00

**Glass beads**

For gentle cleaning and compacting of sensitive surfaces of precious and non-precious alloys. For shining and smoothing out of crowns, bridges and CrCo techniques.

Silicosis-free, made in Germany

<b>656020 C</b>	<b>Glass beads 50 µm – 8 kg</b>	74,50
<b>656021 C</b>	<b>Glass beads 50 µm – 15 kg</b>	105,00
<b>656019 C</b>	<b>Glass beads 125 µm – 8 kg</b>	74,50
<b>656022 C</b>	<b>Glass beads 125 µm – 15 kg</b>	105,00



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15–17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-0  
[export@wassermann-dental.com](mailto:export@wassermann-dental.com) · [www.wassermann.hamburg](http://www.wassermann.hamburg)



[/WassermannDentalmaschinen](#)